User Manual for EFFICACY BOOSTER

HD-272

VAGHEGGI
PHYTocosMETICI

Index

IT - Manuale Utente di EFFICACY BOOSTER.....	2
EN - User Manual for EFFICACY BOOSTER.....	8
FR - Manuel d'utilisation pour EFFICACY BOOSTER.....	13
DE - Benutzerhandbuch für EFFICACY BOOSTER.....	19
ES - Manual de Usuario de EFFICACY BOOSTER.....	26
HU - Használati útmutató az EFFICACY BOOSTER készülékhez.....	33
SK - Používateľská príručka pre EFFICACY BOOSTER.....	39
PT - Manual do Utilizador para EFFICACY BOOSTER.....	46
NL – Gebruikershandleiding van EFFICACY BOOSTER.....	54
PL - Instrukcja użytkownika EFFICACY BOOSTER.....	62
HR - Korisnički priručnik za uređaj EFFICACY BOOSTER.....	69
SI - Uporabniški priročnik za POSPEŠEVALEC UČINKOVITOSTI.....	77
EL - Εγχειρίδιο Χρήστη του EFFICACY BOOSTER.....	85

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EN - User Manual for EFFICACY BOOSTER

HD-272

Distributor: VAGHEGGI SPA

Address: VIA F. PIGAFETTA, 6 36024 NANTO (VI)

Tel.: +39-0444-639600

1. Recommendation of Use

Please read the entire user manual carefully before using the EFFICACY BOOSTER. This will ensure safe and effective use of the device.

1.1 Warnings

1. This product is used for skin care. Do not look directly at the light source for extended periods, as it may cause eye harm.
2. For individuals under the age of 18, please do not use this product.
3. Do not use this product if you are pregnant, preparing for pregnancy, or breastfeeding.
4. Do not use this device on infants.
5. If you experience discomfort or have a concern about its use, discontinue use immediately and consult your doctor.
6. Do not use the device if you notice abnormal operation, visible damage, or if any components appear worn or frayed.
7. Do not use the device in case of metal or plastic implants, surgical intervention areas, ongoing medical conditions, skin lesions, presence of pacemakers, epilepsy, or tattoos.
8. Keep the device away from heat sources, dust, water splashes, humid environments (e.g., bathrooms), high temperatures, and direct sunlight.
9. Use of the device is not recommended in the presence of one of the following conditions:
 - a flammable anaesthetic mixture with air, oxygen or with nitrous oxide.
 - while the user is transit (e.g. during vehicle movement)
 - near magnetic resonance imaging (MRI) equipment.
 - be used adjacent to or stacked with other electronic equipment (unless verified for compatibility)

1.2 Contraindications

1. Do not use this product if you are pregnant, preparing for pregnancy, or breastfeeding.
2. If you have photosensitive skin, please do not use this product.
3. Do not use this device to treat raised, crusted, or bleeding skin.
4. If you have experienced sunburn within the past 2 weeks, do not use this device.
5. Do not use this product within one month of any facial fillers, plastic surgery or other infusions. After one month, this device may be safely used on those treatment areas.
6. Do not use if you are taking drugs known to cause photosensitive reactions, as this may result in skin discoloration, discomfort, or pain.
7. Do not use if you have undergone surgery within the past three months.
8. Do not use if you have undergone radiation therapy or chemotherapy within the past three months or have areas of potential skin malignancies.

1.3 Adverse Reactions

1. This product has no side effects under normal use. However, if you are sensitive to light, please discontinue use immediately.
2. If you experience any discomfort, skin discoloration, skin damage, or any other concerns develop, discontinue use immediately and consult your physician.

1.4 Precautions

1. The device is designed for aesthetic purposes only and for individuals in good health. .If you have any underlying health conditions, consult your physician before using the device.
2. If you experience any discomfort during use, discontinue use immediately.
3. If you suspect an adverse reaction from using the device, consult your doctor immediately.
4. The use of heat- or cold-producing devices, such as ice packs, may impair the performance of device and increase the risk of injury to the patient.
5. Do not scratch or damage the device.
6. Do not attempt to open or repair the device by yourself, as this may cause damage, malfunction, electrical shock, fire or injury. There are no user-serviceable parts inside.

Note:

Disassembling the device will void the warranty.

7. The device requires special precautions regarding electromagnetic compatibility (EMC).
Do not use the device near sources of electromagnetic interference(EMI) such as CT scanners, RFID systems, security screening equipment (e.g. metal detector, EAS), shortwave radios or microwave equipment.
Note: Some sources of Radio Frequency emitters may not be visible to the user. If you experience unexpected behavior of the product, then it is possible some external electromagnetic interference is responsible. You should remove this source and restart the product to restore full operation.
8. Do not store the device in direct sunlight or on hot surfaces. Exposure to excessive heat may damage the device. You should operate, transport and store the device as specified in this manual. Failure to use and maintain the device according to the instructions outlined in this manual will void the warranty.
Risk of strangulation. Keep the charging cable away from children and do not wrap it around the neck.

2 Product Description

The EFFICACY BOOSTER is a hand-held, rechargeable device powered by a Lithium-Ion battery, designed for single-patient use. It offers four integrated functions-red light, micro current, heating and vibrating massage-to enhance facial wellness and visibly helps reduce wrinkles.

The device features a single button to turn it on or off. This device to be used with creams from our product line.

The device automatically shuts off after the recommended 3-minute treatment time.

2.1 Specification

Trade Name		EFFICACY BOOSTER
Model		HD-272
Intended Use		The EFFICACY BOOSTER (Model HD-272) is intended for the treatment of wrinkles for over-the-counter cosmetic use
Design		Hand-held device
Main Unit Weight		90±10 g
Dimensions		Main unit: 67.13*37.4mm±0.2mm Accessory: 3*11.55 mm± 0.2mm
Housing Materials of main unit		PC+PP+SUS304
Power Source		500mAh, 3.7V, lithium battery
Software/Firmware/Mic reprocessor Control?		Yes
Automatic Shut Off		Yes
User Override Control		Yes
Time Range (minutes)		3 minutes per treatment
Wavelengths		Red: 630 ± 10nm
Power density		20~45mW/cm ² ± 10%
Heating Temperature		Max 42°C
Micro current	Wave Form	Pulsed Biphasic Pulsed Biphasic
	Wave Shape	Modulated Square
	Maximum Output Voltage	20VDC

Vagheggi Spa

	Duty Cycle	50%
	Output Tolerance	±10%
Environment for Operating	Temperature: 0 °C~ 50°C/ 32°F~122°F Humidity:30-95% RH Atmospheric Pressure:700h Pa-1060hPa	
Environment for Storage	Temperature: -5 °C~65°C/ 23°F~149°F Humidity:20-95% RH Atmospheric Pressure:700hPa a-1060hPa	
Vibration	Built in	
Charging supply	5v=1A	
Charging time	60±10 minutes	

2.2 Accessories List

The components and accessories included in each package are listed below:

Name	Quantity
Main Unit	1
User Manual	1

2.3 Main Device

ON-Off button

USB-C port

Figure 3-1. Main Device

3 Using Instructions

3.1 First Use and Preparation

1. Remove your device from the package.
2. Connect the USB -C of the charging cable to the device. Plug the USB-A of the charging cable into a suitable USB power adapter (5V DC, 1A output) compliant with regulatory standards.
3. Connect the power adapter to an properly installed electrical outlet (voltage range: 100-240Vac, frequency: 50Hz/60Hz).

Note:

- a. Charge your device for 40 minutes before first use. While charging, the indicator light flashes. It will light remain on when fully charged.
- b. We recommend a full charge cycle of 60 ± 10 minutes if the battery power is drained completely.
- c. If not in use, we recommend charging the device once every three months to maintain battery health.
- d. The device cannot be used while being charged.
- e. The adoption of USB Type-C aims to improve consumer convenience, ensure compatibility with common charging equipment and support environmental sustainability by reducing electronic waste.

3.2 Usage Instructions For Therapy Treatment

3.2.1 Red LED lights and outputs

The device emits red light (wavelength: $630\text{nm} \pm 10\text{nm}$), gentle microcurrent, and heating for the treatment of facial wrinkles.

3.2.2 Treatment operation

1. Clean the area you need to treat before using the device and ensure that the treatment site does not have an open wound.
2. Twist the device to remove the cap.
3. Place the device in the treatment area.
4. Press and hold the power button to turn on the device. Once activated, the red lights will illuminate, the device will begin to vibrate and gradually warm up, and the gentle microcurrent will also be applied to the skin.
5. If your skin feels hot or uncomfortable during the treatment, discontinue use immediately! You might be sensitive to light.
6. The device will automatically shut off after 3 minutes. If you want to stop a treatment early, press the power button again to turn off your device.

4 Maintenance

It is recommended to clean the device after each use.

1. Turn off the device by pressing and holding the power button.
2. Your device requires no maintenance beyond basic cleaning. To clean the device, unplug it and wipe the surface with a damp cloth using water or 70% alcohol. Do not use harsh cleaners or bleach, as they may damage the device. Ensure the device is completely dry before each use.
3. If the cord or device becomes damaged or fails to operate properly, contact the customer support. This device must only be opened, serviced or repaired by VAGHEGGI SPA. All devices returned to the

manufacturer for repair, Whether under warranty or out-of-warranty- must be handled by VAGHEGGI SPA only. Do not attempt to open or repair the device by yourself.

4.1 Trouble Shooting

In the event that the device fails to perform as intended, refer to the following table to identify and resolve potential problems with the device and its setup.

Problem	Solution
Device turns off automatically	Battery
Device cannot turn on	Check and charge. Make sure that the device is in good contact with the charging cable. Make sure the battery is charged.
Device shuts off in mid-treatment or only after few treatments	Check and recharge batteries.
The surface is cracked more than 2mm or shows signs of safety risk due to cleaning	Discontinue use and dispose of the device following the method described in Section 5.3.

4.2 Disposal of Electrical and Electronic Equipment Waste (E-waste)

This device should not be treated as household waste, but rather be brought to the appropriate collection point for recycling of electrical and electronic equipment. The device contains a lithium-ion battery. By ensuring this device is disposed of correctly, you will help prevent the potential negative consequences for the environment and human health which could be caused by inappropriate waste handling of the device. The recycling of materials will also help conserve natural resources.

For more information about the recycling of your device, please contact your local household waste disposal service or your place of purchase.

5

Warranty

The product is guaranteed for a period of 3 years from the date of purchase against manufacturing or material defects.

Claims under warranty must be submitted within the warranty period. Vagheggi, at its sole discretion, will replace the device free of charge. The warranty covers components essential to the correct functioning of the device. It does not cover deterioration due to normal wear and tear, accidents, misuse, or negligence. Any attempt to disassemble the device (or its components) will void the warranty. Warranty service requests must be supported by sufficient proof showing that the request falls within the warranty period.

Therefore, please keep the purchase receipt together with these warranty conditions. In case of malfunction, contact us at info@vagheggi.com with photos and the necessary information

5.1 Outside Warranty Period

Any damage caused by non-original packaging may render the warranty null and void. Any services to the device shall be performed only by a service technician certified by VAGHEGGI SPA.

Vagheggi Spa

6 Contact

	IP Classification: IP22 ; Protected against solid foreign objects $\geq 12,5$ mm and against vertically falling water drops when enclosure tilted up to 15°
	Indicates the device is CE marked and authorized to sell in Europe. (Lower-Voltage-Directive)
	The Waste of Electrical and Electronic Equipment (WEEE). Indicates the device may not be discarded as household waste (2012/19/EU)
	Indicates that the device complies with the Restriction of Hazardous Substances Directive (2011/65/EU) and does not contain substances such as lead, mercury, cadmium, or other restricted hazardous materials above regulated thresholds.
	Indicates that the device complies with the REACH Regulation (EC 1907/2006) concerning the Registration, Evaluation, Authorization and Restriction of Chemicals, and does not contain Substances of Very High Concern (SVHC) above the legal limits.

Information:

Importer: VAGHEGGI SPA

Address: VIA F. PIGAFETTA,6 36024 NANTO (VI)

Tel: +39-0444-639600

Authorized European representative: Light Tree Ventures Holding B.V.
 Laan van Ypenburg 108, 2497 GC, The Hague, The Netherlands
regulatory@lighttreeventures.com

Manufacturer: Shenzhen Nuon Medical Equipment Co., Ltd.

Address: 1st Floor-3rd Floor, No. 27-2,Xintang Road,Xintian Community, Fuhai Street, Baoan District,Shenzhen, Guangdong, 518103, China

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